



EliZyme™ HS HIFI MIX

Intended use:

For Research Use Only. Not for use in diagnostic procedures.

Storage:

Upon arrival store components at -20 °C. Avoid prolonged exposure to light. When stored under these conditions, the kit will retain full activity until the expiration date indicated on the kit label. Reagents may be stored at 4 °C up to 1 month.

Product description

EliZyme™ HS HIFI MIX has been engineered to exhibit a greater affinity for DNA, eliminating the need for accessory proteins or DNA binding domains. This enzyme possesses a natural high processivity, leading to notable enhancements in yield, speed, and sensitivity. It demonstrates improved capabilities in amplifying long DNA fragments, as well as targets with high GC- or AT-richness. To prevent non-specific amplification during the setup of the reaction, enhance sensitivity, and improve reaction efficiency, the enzyme is combined with a proprietary antibody that deactivates it until the initial denaturation step. EliZyme™ HS HIFI MIX is supplied in a convenient 2X concentrated format. The 2X EliZyme™ HS HIFI MIX comprises HotStart DNA Polymerase, dNTPs, MgCl₂, and stabilizers. The EliZyme™ HS HIFI MIX is specifically designed for routine, high-fidelity PCR of various targets and fragment sizes. It offers error rates that are approximately 100 times lower than those of the standard Taq DNA polymerase.

EliZyme™ HS HIFI polymerase exhibits 5'-3' polymerase activity and 3'-5' exonuclease (proofreading) activity, but lacks 5'-3' exonuclease activity. The strong 3'-5' exonuclease activity contributes to extremely accurate DNA amplification. The error rate of EliZyme™ HS HIFI polymerase is 1 error per 3.6 x 10⁶ nucleotides incorporated. DNA fragments produced using the EliZyme™ HS HIFI MIX are suitable for routine downstream analysis and applications, such as restriction enzyme digestion, cloning, and sequencing. PCR products have blunt ends.

Content

| | Ref. No. | Content | Size |
|----------------------|----------|-----------|-----------|
| EliZyme™ HS HIFI MIX | EZ2501 | 1x1.25 ml | 100 rxns |
| | EZ2505 | 1x6.25 ml | 500 rxns |
| | EZ2510 | 2x6.25 ml | 1000 rxns |



Primers

Primers should have a predicted melting temperature of around 65 °C. Primers should be designed to eliminate the possibility of primer-dimer formation and non-specific amplification. The final primer concentration in the reaction should be between 0.2 µM and 0.6 µM.

Reaction setup

After thawing, briefly vortex the mix and shortly spin.

| Reagent | 25 µl reaction | Final conc. |
|-------------------------|------------------------------------|-------------|
| 2x EliZyme™ HS HIFI MIX | 12.5 µl | 1x |
| Forward primer (10 µM) | 0.75 µl | 300 nM |
| Reverse primer (10 µM) | 0.75 µl | 300 nM |
| Template DNA | < 100 ng genomic DNA, < 10 ng cDNA | Variable |
| PCR grade water | Up to 25 µl | |

PCR cycling profile

| Step | Temperature | Time | Cycles |
|----------------------|-------------|-------------|----------|
| Initial denaturation | 95 °C | 1 – 3 min | 1 |
| Denaturation | 98 °C | 20 s | |
| Annealing | 60 – 75 °C | 15 s | 25 – 35* |
| Extension | 72 °C | 15 – 60 s** | |
| Final extension | 72 °C | 1 min/kb*** | 1 |

*For highest fidelity is possible to use less than 25 cycles.

**For fragments longer than 1 kb use up to 60 s/kb.

***Optional.

Manufacturer:

ELISABETH PHARMACON, spol. s r. o.

Rokycanova 4437/5, Brno-Židenice 615 00

info@elisabeth.cz | www.elisabeth.cz | tel.: +420 542 213 851



Catalog number



Batch code



Use by (last day of month)



Upper limit of temperature



Manufacturer



Contains sufficient "N" tests