

# Bead Ruptor 12

## BEAD MILL HOMOGENIZER



Compelling performance  
at a compelling price



# Bead Ruptor 12



ORDER NO. 19-050A

**Ideally suited for nucleic acid, drug, protein, enzyme, and virus extractions.**

## Specifications

**Processing Range:** 25  $\mu$ L to 6 mL

**Power Rating:** 115-230 VAC/50-60 Hz

**Dimensions:** 13" W x 15.5" D x 13" H

**Weight:** 29 kg (65 lbs.)

**Compliance:** CE

**Run Time:** 1 s - 9:59 m

**Performance Range:** 0.8 m/s - 6 m/s  
in increments of 0.15 m/s

**Cycles:** 1-10

**Programs:** 99


## Includes Tube Carriage for 2 mL and 7 mL tubes

- Holds 12 x 2 mL or 0.5 mL tubes and/or 4 x 7 mL tubes
- Includes semi-permanent ring



## Features

- Optimized tube motion with speeds up to 6 m/s
- Process up to 16 samples per cycle
- Convenient front loading design
- Locking safety lid
- Reproducible extraction of proteins, nucleic acids, and small molecules



935 Cobb Place Blvd.  
Kennesaw, GA 30144  
800.776.4431 • 770.421.0058  
www.omni-inc.com

REV.D