



✓ DESCRIPTION

EliGene® Urine Isolation Kit is intended for the isolation of DNA from clinical material.

EliGene® Sperm Isolation Supplement is intended for DNA isolation from the sperm.

Many problems like inhibitions, low sensitivity etc. may affect isolation of DNA from clinical samples. **EliGene® Urine Isolation Kit** is designed to avoid all problems with PCR inhibitors and low yields of DNA. The total time of isolation process is about 1 – 2 hour including incubation. Kit contains all components necessary for the DNA isolation including microtubes. No addition of ethanol or other chemicals is needed.


✓ PERFORMANCE STUDY

Within the frame of testing, the functional characteristics of EliGene® Chlamydia trachomatis LC kit overall 120 clinical specimens of urine were analysed. Before the detection by EliGene® Chlamydia trachomatis LC, DNA from all samples was isolated by the EliGene® Urine Isolation Kit. No inhibition was observed. The effectiveness of this isolation procedure was 100%. It means that from total amount of isolated specimens were 100% of specimens without inhibition.

✓ ADVANTAGES

- ✘ Rapid process of DNA isolation from urine in less than 20 minutes
- ✘ High quality isolated DNA suitable for any downstream application (PCR, qPCR, restriction digest)
- ✘ Complete package - all required components included in the kit

✓ EASY PROTOCOL



CONTACT US

ELISABETH PHARMACON
ROKYCANOVA 4437/5, BRNO-ZIDENICE 615 00
CZECH REPUBLIC | +420 542 213 851
orders@elisabeth.cz | www.elisabeth.cz

