



✓ DESCRIPTION

EliGene[®] Celiac DQ PLUS RT is intended for the screening of alleles

- HLA-DQ2.5 (DQA1*05/DQB1*02),
- HLA-DQ2.2 (DQA1*02/DQB1*02),
- HLA-DQ8 (DQA1*03/DQB1*03:02)
- HLA-DRB1*04 from isolated DNA.

Specimens: Buccal swabs, blood
 Detection technology: TaqMan LNA

✓ PERFORMANCE STUDY

EliGene[®] Celiac DQ PLUS RT kit specificity was tested on 50 samples of human DNA with genotypes determined by CE-IVD kits for the detection of analysed alleles. Totally 50 samples were determined right by EliGene[®] Celiac DQ PLUS RT

Clinical specificity: 100%.

Analytical sensitivity: from 5 ng DNA molecules.


✓ ORDER DETAILS

Cat. No.: 90085-RT

50 reactions

✓ FEATURES

- Detection of all HLA alleles associated with Coeliac Disease
- Discrimination of DQ2.2, DQ2.5, DQ8 and DRB1*04 alleles
- Multiplex Real-Time PCR solution with results in 60 min
- In accordance with the latest recommendations of gastroenterological experts
- Including internal and positive control



✓ COMPATIBLE DEVICES

QuantStudio 5 (ThermoFisher Scientific), RotorGene Q (Qiagen), CFX96 Touch (Bio-Rad), MIC Cycler (Bio Molecular Systems), etc.

CONTACT

ELISABETH PHARMACON
 ROKYCANOVA 4437/5, BRNO-ŽIDENICE
 CZECH REPUBLIC | +420 542 213 851
www.elisabeth.cz

