



✓ DESCRIPTION

EliGene® Coeliac 3.0 RT is intended for the screening of alleles

- HLA-DQ2.5 (DQA1*05/DQB1*02),
 - HLA-DQ2.2 (DQA1*02/DQB1*02),
 - HLA-DQ8 (DQA1*03/DQB1*03:02)
- from isolated DNA.

Specimens: Buccal swabs, blood

Detection technology: TaqMan LNA

✓ PERFORMANCE STUDY

EliGene® Coeliac 3.0 RT kit specificity was tested on 50 samples of human DNA with genotypes determined by DNA sequencing. Totally 50 samples were right determined by EliGene® Coeliac 3.0 RT.

The clinical specificity of EliGene® Coeliac 3.0 RT kit is 100%.

Analytical sensitivity: from 1 ng DNA molecules.

Specificity: 100 %


✓ ORDER DETAILS

Cat. No.: 90071-RT

50 reactions

✓ ADVANTAGES

- Detection of all HLA alleles associated with Coeliac Disease
- Discrimination of DQ2 suballeles: DQ 2.2 and DQ 2.5
- Real-Time PCR
- Results in 60 min.
- In accordance with the latest recommendations of gastroenterological experts
- Including internal and positive control



✓ COMPATIBLE DEVICES

ABI 7500 FAST (Applied Biosystems); LightCycler®480 (Roche), RotorGene Q (Qiagen); Quant Studio 3 and 5 (Thermofisher Scientific); CFX96 Touch Real -Time PCR Detection System (Bio -Rad)

CONTACT

ELISABETH PHARMACON
 ROKYCANOVA 4437/5, BRNO-ŽIDENICE
 CZECH REPUBLIC | +420 542 213 851
www.elisabeth.cz

